

12:04

Your Nightscout site and API

https://testowycgm.herokuapp.com

test

SAVE

ENABLE AUTO STOP

AUTO STOP AT BG VALUE: 120MG/DL

BG SOURCE: PUMP / NS

ENTER PHONE NUMBER

SEND LOGS

OK! https://testowycgm.herokuapp.com

NIGHTSCOUT

BLUETOOTH

CONFIG

12:44

Bluetooth and pump operations

BLUETOOTH SCAN FOR MED-LINK

254

MED-LINK ,90:E2:02:8D:62:21

RUN AND ACCEPT THE TERMS:

MAKE RESUME PUMP

Aug 24 2020 12:44  
WARNING!  
PUMP SUSPENDED BY AUTO STOP FUNCTION  
  
24-08-2020 12:43 58%  
BG: 98 24-08-20 12:39  
Last bolus: 0.6u 24-08-20 12:05  
Square bolus: 0.0u delivered: 0.000u  
Square bolus time: 0h:00m / 0h:00m  
ISIG: 23.80nA  
Calibration factor: 4.485  
Next calibration time: 0:09  
Sensor uptime: 34min  
BG target: 75-155  
Pump battery voltage: 1.48V  
Reservoir: 200.77u  
Reservoir: 200.77u

NIGHTSCOUT

BLUETOOTH

CONFIG

12:47

Bluetooth and pump operations

BLUETOOTH SCAN FOR MED-LINK

253

MED-LINK ,90:E2:02:8D:62:21

RUN AND ACCEPT THE TERMS:

MAKE SUSPEND PUMP

Aug 24 2020 12:47  
WARNING!  
PUMP RESUMED BY AUTO START FUNCTION  
  
24-08-2020 12:46 57%  
BG: 106 24-08-20 12:44  
Last bolus: 0.6u 24-08-20 12:05  
Square bolus: 0.0u delivered: 0.000u  
Square bolus time: 0h:00m / 0h:00m  
ISIG: 25.70nA  
Calibration factor: 4.485  
Next calibration time: 0:09  
Sensor uptime: 37min  
BG target: 75-155  
Pump battery voltage: 1.48V  
Reservoir: 200.77u  
Reservoir: 200.77u

NIGHTSCOUT

BLUETOOTH

CONFIG